



PRODUCT CODE: RE-SURMA-3PLY

PRODUCT DESCRIPTION: Disposable 3-PLY surgical mask

MATERIALS

NOSECLIP	Aluminium strip, glue bonded
MASK/RESPIRATOR	3 layers, top and bottom are non-woven 100% polyester, middle layer is meltblown nonwoven polyester
EAR STRAPS	Spandex elastic tape, glue bonded to mask

COMPLIANCY

EN/ISO REQUIREMENT	<ul style="list-style-type: none"> • YY0469-2011 • GB1579-2002 ANNEX B • EN14683
SPECIFIC TESTS	<ul style="list-style-type: none"> • Bacterial Filtration Efficiency (BFE) – all tests > 99.5 • Packed in a sterile environment

SPECIAL FEATURES

- Sterile surgical mask, 3-PLY. Tested to YY0469-2001 or EN14683 standards.
- Pre-stretched, latex-free head straps.

WARNINGS

This product does not supply oxygen. Ensure that you use it in well ventilated areas containing sufficient oxygen.

This is a single use respirator. Discard it after a single shift, or after 8 hours of use.

Leave work area immediately if breathing becomes difficult, or dizziness occurs.

Consult your Safety Officer or an approved sales rep when implementing a respiratory protection programme.

OTHER INFORMATION

SIZING	20 masks in a single packet
PACKAGING	2000 masks in a carton
STORAGE	Respirators must be placed in a sealed packet or container, and stored in an environment free of dust, chemicals, moisture, excessive heat and direct sunlight. Failure to do so may affect the performance and shelf life of these products.